



CUSTOMER CODE 7193/TQC163-13

Specific Proteins 05/2016

Sample SR0216

TQSR0030-May-2016

C3 (mg/dl)		Instrument					
		Reagent					
		Method: Turbidimetry- Roche (own group)					
		Observations					
Distribution of the results		Results history:					
<p>Own lab result: 55,00 mg/dl</p>		Month	Sample code	X (Mean)	Results	Bias (Diff%)	Target-limit
		5/2016	SR0216	55,46	55,00	-0,82	
		3/2016	SR0116	1,06	1,14	8,06	
		11/2014	SR0414	105,11	110,00	4,66	
		9/2014	SR0314	105,59	110,00	4,18	
		5/2014	SR0214	54,00	54,00	0,00	
		3/2014	SR0114	111,00	111,00	0,00	
		11/2013	SR04	55,00	55,00	0,00	
		9/2013	SR03	102,74	99,00	-3,64	
		5/2013	02	0,91	0,88	-3,20	
		3/2013	01	0,62	0,62	0,00	
			-	-	-	-	
	-	-	-	-			
						-30% x +30%	
Results	x (mean)	SD (standard deviation)	CV%	n (no. of results)	x±30%(target-limit)	Z score	
Own group	55,46	1,39	2,50	7	38,82 - 72,09	-0,33	
All methods	64,48	32,40	50,25	69	45,14 - 83,83	-0,29	