



**Laboratory Disposables**




**Medical Disposables**

# Laboratory & Medical Disposables



[www.todylaboratories.com](http://www.todylaboratories.com)



**CONTAINERS FOR SPECIMEN**

**PP, transparent polypropylene, with polyethylene red screw cap**

urine container, graduated at every 5 ml till 30 ml, frosted label	30 ml	C22LD1005
urine container, individually wrapped, graduated at every 10 ml till 30 ml, total vol. 40 ml, frosted label	30 ml	C51LD1032
urine container, sterile, individually wrapped, graduated at every 10 ml till 30 ml, total vol. 40 ml, frosted label	30 ml	C51LD1033
urine container, sterile, individually wrapped, graduated at every 5 ml till 30 ml, frosted label	30 ml	C22LD1006
faeces container, sterile, individually wrapped, graduated at every 5 ml till 30 ml, frosted label	30 ml	C22LD1003



**PP, transparent polypropylene**  
**PP, black polypropylene**

faeces container, with attached spoon, neutral	20ml	C51LD1030
faeces container, with attached spoon, black, recommended for photosensitive samples	20ml	C51LD1031


**PP, transparent polypropylene, with polyethylene blue screw cap**  
**PP, transparent polypropylene, with polyethylene red screw cap**


urine container, sterile, individually wrapped, graduated at every 10ml till 60 ml, with label	60 ml	C51LD1021
urine container, sterile, individually wrapped, graduated at every 5 ml till 60 ml, with frosted label	60 ml	C22LD1007
urine container, sterile, individually wrapped, with frosted label, graduated every 10 ml till 100 ml	120 ml	C22LD1004

**INOCULATING LOOPS**

**PS, polystyrene, white**  
**PS, polystyrene, light green**  
**PS, polystyrene, dark blue**  
**PS, polystyrene, light blue**

soft, sterilized by Gamma Rays	1µl	C51LD1002
soft, sterilized by Gamma Rays	1µl	C22LD1001
soft, sterilized by Gamma Rays	10µl	C51LD1003
soft, sterilized by Gamma Rays	10µl	C22LD1002





**TOURNIQUET**  
Terylene elastic belt

ABS auto lock buckle, blue

450 x 250mm

**code**  
C51LD1034

**MICROSCOPE SLIDES**  
Glass

frosted 1 end 1 side, 90° corners  
cut edges, plain, 90° corners  
frosted 1 end 2 sides, 45° clipped corners

75x 25x 1,1mm  
75x 25x 1,1mm  
75x 25x 1,1mm

C51LD1005  
C51LD1006  
C51LD1007



**MICROSCOPE COVERGLASS**  
Glass

thickness 0,13-0,16 mm  
thickness 0,13-0,16 mm  
thickness 0,13-0,16 mm  
thickness 0,13-0,16 mm  
thickness 0,13-0,16 mm

18 x 18 mm  
20 x 20 mm  
22 x 22 mm  
24 x 24 mm  
24 x 50 mm

C51LD1011  
C51LD1009  
C51LD1008  
C51LD1010  
C51LD1012




**PASTEUR PIPETTE**  
LDPE, transparent polyethylene

total capacity 7 ml, graduated at every 0,5 ml till 3 ml, 3 ml  
length: 155 mm; not sterile

C51LD1004





**PASTEUR PIPETTE**

Transparent polyethylene

total capacity 7 ml, graduated at every 0,5 ml till 3 ml, 3 ml  
length: 155 mm; not sterile

**code**  
C22LD1008

**STERILE, LDPE, transp. polyeth.,  
graduated, individually flow-pack**

graduated, length: 150 mm; sterile 5 ml

C51LD1027

**PETRI DISHES**

PS, polystyrene,  
high transparency

triple vents, not sterile, produced in cleanroom ISO 6 Ø 90 mm  
h 16,2 mm

C22LD1010

**STERILE, PS, polystyrene,  
high transparency**

sterile, suitable with the most automated Petri dish  
filling instruments Ø 90 mm  
h 15 mm

C51LD1028



**CYLINDRICAL TEST TUBES**

PP, transp. polypropylene

12 x 75 mm 5 ml

C51LD1013

13 x 75 mm, suitable with stopper C51LD1022 5 ml

C51LD1014

**STERILE, PP, transparent  
polypropylene, suitable with  
stopper C51LD1024**

16 x 100 mm 10 ml

C51LD1029

**PP, transp. polypropylene,  
suitable with stopper C51LD1024**

16 x 100 mm 10 ml

C51LD1025



**MULTIWELL**

96 wells

Multiwell (microplate), U-shape bottom, without lid, 96  
PS material


C51LD1035

**CONICAL CENTRIFUGE TUBES**

**STERILE, PP, transparent  
polypropylene, with screw  
gray lid, graduated**

17 x 120 mm 15 ml

C51LD1026



**MEDICAL SWABS**

**Stick + cotton tip**

		code
wooden stick with cotton tip length, individually wrapped, sterile	150 mm	26075-500
wooden stick with cotton tip in polypropylene test tube, sterile	12 x 150 mm	26060-100
wooden stick with cotton tip in polypropylene test tube, sterile, sterilised with Gamma rays	13 x 150 mm	2120-0010
wooden stick with cotton tip in polypropylene test tube, sterile, sterilised with E.O. rays	13 x 150 mm	2120-0009
individually packed, sterilised with E.O. rays	2.2 x 150 mm	2122-0002
individually packed, non-sterile	2.2 x 150 mm	2122-0007
individually packed, sterile	2.2 x 150 mm	2122-0009
in polypropylene test tube with Cary Blair transport medium, individual packed, sterile	13 x 150 mm	2120-0028

**PS, polystyrene stick + viscose tip**

**PS, polistiren + viscose tip with Cary Blair transport medium**

**Plastic stick with viscose tip, in polypropylene test tube with Amies transport medium**

plastic stick with rayon tip in polypropylene test tube, sterile, used for transportation of various pathogens like Neisseria sp., Haemophilus sp., Corynebacterium diphteriae, Trichomonas vaginalis, Streptococcus pyogenes, Streptococcus pneumoniae and Salmonella	13 x 150 mm	2120-0026
	12 x 150 mm	26068-100



**CONICAL MICROTUBES**

**PP, transp. polypropylene**



10.6 x 39.5 mm, non-sterile graduated at: 0,5-1-1,5 ml, with flat cap and frosted label, one-handed usage, rezistant at centrifugation 15.000G and temperatures from -86°C to +121°C	1,5 ml	C51LD1015
graduated at: 0,5-1-1,5 ml, with flat cap and frosted label, one-handed usage, rezistant at centrifugation 15.000G and temperatures from -70°C to +120°C	1,5 ml	C22LD1016

**MICROTUBES ROUND BOTTOM**

**PP, transp. polypropylene**

10.3 x 40.5 mm, autoclavable, with flat cap and frosted label, one-handed usage, rezistant at centrifugation 10.000G and temperatures from -86°C to +121°C	2 ml	C51LD1016
--	------	-----------

**TIPS FOR PIPETTES**

**PP, polypropylene**



0-200 µl, Gilson type, bevelled	yellow	C51LD1018
	yellow	C22LD1019
0-200 µl, Eppendorf, Brand, Socorex type, with collar	yellow	C51LD1019
	yellow	C22LD1018
0,5-10 µl, Eppendorf type	neutral	C51LD1017
201-1000 µl, Eppendorf type	light blue	C51LD1020
1000 µl, Gilson type, without collar	light blue	C51LD1023
201-1.000 µl, Eppendorf, Brand, Socorex Gilson type, w/collar	light blue	C22LD1017



**RIGID FROZEN PACK -21°C  
COLLAPSIBLE ISOTHERMAL BOX  
- LIPBOX**

900 ml, 180x143x36mm, 753gr  
2,7 litres, 205x118x110mm, thickness 4mm, FXPS  
12 litres, 300x200x200mm, thickness 4,5mm, FXPS  
18 litres, 320x235x240mm, thickness 5mm, FXPS  
30 litres, 400x300x250mm, thickness 4,5mm, FXPS  
54 litres, 420x315x410mm, thickness 4,5mm, FXPS

1-10137/RFP900  
1-RM35/FXPS-3-02  
1-DISF12/FXPS-12-01  
1-NOE02/FXPS-18-01  
1-DISF30/FXPS-30-01  
1-IV01/FXPS-54-01  
1-DISF12-DBL/FXPS-12-DBL

**code**

**COLLAPSIBLE ISOTHERMAL BOX,  
DOUBLE WALL - LIPBOX**


12 litres, 300x200x200mm,  
thickness 9mm = 2x4,5mm, FXPS  
18 litres, 320x235x240mm,  
thickness 9mm = 2x4,5mm, FXPS  
30 litres, 400x300x250mm,  
thickness 9mm = 2x4,5mm, FXPS  
32,9 litres, 415x315x250mm,  
thickness 9mm = 2x4,5mm, FXPS

1-NOE02-plus-DBL/FXPS-18-DBL  
1-DISF30-DBL/FXPS-30-DBL  
1-MIL-02-NII/FXPS-31-01DBL

**PU PACKAGING**

15 litres, external: 390x295x370mm, internal:  
290x195x268mm; 1.339gr; thickness 40mm  
29 litres, external: 471x401x362mm, internal:  
371x301x263mm; 2.000gr; thickness 40mm  
29 litres, external: 510x435x410mm, internal:  
371x295x265mm; 2.580gr; thickness 60mm  
40 litres, external: 510x495x440mm, internal:  
371x360x300mm; 3.540gr; thickness 60mm  
48 litres, external: 440x440x525mm, internal:  
340x340x420mm; 2.540gr; thickness 40mm  
71 litres, external: 568x458x522mm, internal:  
470x360x420mm; 3.520gr; thickness 40mm  
96 litres, external: 705x505x505mm, internal:  
600x400x400mm; 4.000gr; thickness 40mm

CIPU40-15L-CCE-TPS/CIPU-15P  
CIPU40-29L-CCE-TPS/CIPU-29P  
CIPU60-29L-CCE-TPS/CIPU-29N  
CIPU60 40L-CCE-TPS/CIPU-40N  
CIPU40-48L-CCE-TPS/CIPU-48P  
CIPU40-71L-CCE-TPS/CIPU-71P  
CIPU40-96L-CCE-TPS/CIPU-96P



**code**

**REUSABLE ISOTHERMAL XPS - LIPCONTAINER**



4,6 litres, external: 312x279x314mm, internal: 177x145x180mm; 1510gr; thickness 60mm, for frozen products

7,5 litres, external: 384x274x238mm, internal: 291x181x143mm; 2.200gr; thickness 40mm, for refrigerated or Environmental (+25° C) products

13,3 litres, external: 391x291x363mm, internal: 291x181x253mm, 2.830gr; thickness 40mm, for refrigerated or Environmental (+25° C) products

20,8 litres, external: 504x430x324mm, internal: 370x295x190mm; 3.140gr; thickness 60mm, for frozen products

29 litres, external: 504x430x399mm, internal: 371x295x265mm; 3.220gr; thickness 60mm, for frozen products

29,4 litres, external: 463x394x358mm, internal: 371x301x263mm, thickness 40mm, for refrigerated or Environmental (+25° C) products

40 litres, external: 504x495x435mm, internal: 371x360x300mm; 4.040gr; thickness 60mm, for frozen products

74,4 litres, external: 575x487x522mm, internal: 477x377x414mm; thickness 40mm, for refrigerated or Environmental (+25° C) products

96 litres, external: 698x499x510mm, internal: 600x400x400mm; 4400gr; thickness 40mm, for refrigerated or Environmental (+25° C) products

**BAGS WITH GEL TPS THERMAL BARRIER**

220x300x10mm, 310gr, thickness 10mm, vol. 250 cc  
220x300x20mm, 555gr, thickness 20mm, vol. 500 cc

325x453x10mm, 1.150gr, thickness 10mm, vol. 850 cc  
350x300x10mm, 435gr, thickness 10mm, vol. 350 cc

350x300x20mm, 645gr, thickness 20mm, vol. 600 cc  
510x270x10mm, 645gr, thickness 10mm, vol. 650 cc

510x270x20mm, 1.100gr, thickness 20mm, vol. 850 cc  
520x560x10mm, 1487gr, thickness 10mm, vol. 1000 cc

520x560x20mm, 1.674gr, thickness 20mm, vol. 1000 cc  
size S (small), 250x285mm, 250 gr, double wall

size M (medium), 300x350mm, 375 gr  
size B (big), 400x485mm, 500 gr

**BAGS ISOTERM WITH HANDLE**

**REUSABLE ISOTHERMAL - DIAGNOSACH**

5,56 litres, bag: 250x250x290mm; container: 185x165x182mm, 0.95kg, thickness 30mm

12,8 litres, bag: 400x300x300mm; container: 325x185x213mm, 1,55kg, thickness 40mm

22,6 litres, bag: 450x350x350mm; container: 355x255x250mm, 1,98kg, thickness 40mm

42,2 litres, bag: 450x450x430mm; container: 355x355x335mm, 3,4kg, thickness 40mm

1-TB-00N-KF/CIXPS-5N

1-TB-01P-KF/CIXPS-7P

1-TB-02P-KF/CIXPS-13P

1-TB-01N-KF/CIXPS-21N

1-TB-02N-KF/CIXPS-29N

1-TB-03P-KF/CIXPS-29P

1-TB-03N-KF/CIXPS-40N

1-TB-04P-KF/CIXPS-74P

1-TB-05P-KF/CIXPS-96P

1-KF-FLEX-325/TB250S

1-KF-FLEX-325D/TB250D

1-KF-FLEX-860S/TB850S

1-KF-FLEX-720/TB350S

1-KF-FLEX-720D/TB350D

1-KF-FLEX-SPM6S/TB650S

1-KF-FLEX-SPM6D/TB650D

1-KF-FLEX-SPM16S/TB1000S

1-KF-FLEX-SPM16D/TB1000D

1-10524/COLBAG25x28

1-10525/COLBAG30x35

1-10528/COLBAG40x48

1-DIAGNOSACH-425-KF/DSX-5

1-DIAGNOSACH-12-KF/DSX-12

1-DIAGNOSACH-21-KF/DSX-21

1-DIAGNOSACH-40-KF/DSX-40



🏢 **Tody Laboratories Int. SRL**  
🏠 22 Vadul Moldovei, s1, OP 18  
Bucharest, CP 014033, ROMANIA  
📞 ISO 15189, 17043, 9001, 13485  
🌐 [www.todylaboratories.com](http://www.todylaboratories.com)

**tody**laboratories®  
Responsible diagnostic. Since 1993.

[www.todylaboratories.com](http://www.todylaboratories.com)

