



Vacuum Tubes

Collection Needles


Injection Needles

2024

VACUUM TUBES NEEDLES FOR VACUUM TUBES AND SYRINGES



www.todylaboratories.com




CLINICAL CHEMISTRY & IMMUNOLOGY



PRODUCT NAME	DESCRIPTION	SIZE	CODE
Activating gel	13 x 75 mm, plastic, yellow cap	3 ml	C78VT1004
	13 x 100 mm, plastic, yellow cap	5 ml	C78VT1005
	13 x 100 mm, STERILE, plastic, yellow cap	5 ml	C78VT1024
	16 x 100 mm, plastic, yellow cap	8 ml	C78VT1006
Activating gel - PEDIATRIC USE	8 x 45 mm, plastic, yellow cap	0.5 ml	C78VT1042



PRODUCT NAME	DESCRIPTION	SIZE	CODE
Lithium heparin	13 x 75 mm, plastic, green cap	4 ml	C78VT1034
	13 x 100 mm, STERILE, plastic, green cap	6 ml	C78VT1035




PRODUCT NAME	DESCRIPTION	SIZE	CODE
Lithium heparin - PEDIATRIC USE, micro tubes with integrated capillary	8 x 45 mm, plastic, green cap	0.5 ml	C78VT1046



PRODUCT NAME	DESCRIPTION	SIZE	CODE
Siliconated (activating clot)	13 x 75 mm, plastic, red cap	4 ml	C78VT1001
	13 x 100 mm, plastic, red cap	6 ml	C78VT1002
	16 x 100 mm, plastic, red cap	8 ml	C78VT1003
Siliconated (activating clot) - PEDIATRIC USE	8 x 45 mm, plastic, red cap	0.5 ml	C78VT1043
Siliconate (plain) - PEDIATRIC USE	8 x 45 mm, plastic, red cap	0.5 ml	C78VT1047

GLUCOSE

PRODUCT NAME	DESCRIPTION	SIZE	CODE
Potassium oxalate/Sodium fluoride - PEDIATRIC USE, micro tubes with integrated capillary	8 x 45 mm, plastic, grey cap	0,5 ml	C78VT1048



COAGULATION


PRODUCT NAME	DESCRIPTION	SIZE	CODE
Sodium citrate 3.2%	13 x 75 mm, plastic, blue cap	1.8 ml	C78VT1013
	13 x 75 mm, plastic, blue cap	2.7 ml	C78VT1014
	13 x 75 mm, plastic, blue cap	4.5 ml	C78VT1015



HEMATOLOGY



PRODUCT NAME	DESCRIPTION	SIZE	CODE
EDTA K2	13 x 75 mm, plastic, purple cap	2 ml	C78VT1007
	13 x 75 mm, plastic, purple cap	3 ml	C78VT1008
	13 x 75 mm, plastic, purple cap	4 ml	C78VT1009



EDTA K2 - PEDIATRIC USE	8 x 45 mm, plastic, purple cap	0.5 ml	C78VT1044
EDTA K3	13 x 75 mm, plastic, purple cap	2 ml	C78VT1010
	13 x 75 mm, plastic, purple cap	3 ml	C78VT1011
	13 x 75 mm, plastic, purple cap	4 ml	C78VT1012
EDTA K3 - PEDIATRIC USE	8 x 45 mm, plastic, purple cap	0.5 ml	C78VT1045

ESR


PRODUCT NAME	DESCRIPTION	SIZE	CODE
ESR , reading on a rack with scale. FREE : at each 10 boxes x 100 tubes bought we grant for free a rack	8 x 120 mm, glass, black cap	1.6 ml	C78VT1016
	8 x 120 mm, glass, black cap	1.8 ml	C78VT1017
	8 x 120 mm, glass, black cap	2.0 ml	C78VT1018
ESR rack with scale, 10 positions	-	-	C66VT1019



OTHERS

PRODUCT NAME	CODE
Holders , for collection needles 18G, 20G, 21G, 22G and needles flashback or visio type 21G, 22G	EL310






NEEDLES for VACUUM TUBES

PRODUCT NAME	DESCRIPTION	SIZE	CODE
Needles , sterile, non-pyrogenic, non-toxic, single use, safe and efficient, latex-free, allow multiple samples to be collected at a single puncture, sharp and smooth edges make painless penetration	1.20 mm x 38 mm, pink	18G	ELB18
	0.90 mm x 90 mm, yellow	20G	ELB20
	0.80 mm x 38 mm, green	21G	ELB21
	0.70 mm x 90 mm, black	22G	ELB22





NEEDLES for SYRINGES

PRODUCT NAME	DESCRIPTION	SIZE	CODE
Needles, sterile, disposable, non-toxic, non-pyrogenic	1.27 mm x 38 mm, pink	18G	ELHN18
	0.90 mm x 38 mm, yellow	20G	ELHN20
	0.81 mm x 38 mm, green	21G	ELHN21
	0.71 mm x 32 mm, black	22G	ELHN22
	0.64 mm x 32 mm, blue	23G	ELHN23
	0.55 mm x 38 mm, purple	24G	ELHN24
	0.50 mm x 38 mm, orange	25G	ELHN25
	0.45 mm x 16 mm, brown	26G	ELHN26





🏢 **Tody Laboratories Int. SRL**
🏠 22 Vadul Moldovei, s1, OP 18
Bucharest, CP 014033, ROMANIA
📞 ISO 15189, 17043, 9001, 13485
🌐 www.todylaboratories.com

todylaboratories®
Diagnostic responsabil. Din 1993.

www.todylaboratories.com

