



Vacuum Tubes

Collection Needles


Injection Needles

2025

**VACUUM TUBES
NEEDLES FOR VACUUM
TUBES AND SYRINGES**



www.todylaboratories.com




CLINICAL CHEMISTRY & IMMUNOLOGY



PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Activating gel	13 x 75 mm, plastic, yellow cap	3 ml	100	C75VT1004
	13 x 100 mm, plastic, yellow cap	5 ml	100	C75VT1005
	16 x 100 mm, plastic, yellow cap	8 ml	100	C75VT1006
Activating gel - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, yellow cap	0,5 ml	100	C75VT1042



PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Lithium heparin	13 x 75 mm, plastic, green cap	4 ml	100	C75VT1034
	13 x 100 mm, plastic, green cap	6 ml	100	C75VT1035
Lithium heparin - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, green cap	0,5 ml	100	C75VT1046




PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Siliconated (activating clot)	13 x 75 mm, plastic, red cap	4 ml	100	C75VT1001
	13 x 100 mm, plastic, red cap	6 ml	100	C75VT1002
	16 x 100 mm, plastic, red cap	8 ml	100	C75VT1003
	16 x 100 mm, plastic, red cap	9 ml	100	C75VT1049
	16 x 100 mm, plastic, red cap	10 ml	100	C75VT1050
Siliconated (clot activator) - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, red cap	0,5 ml	100	C75VT1043
Siliconated (plain) - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, red cap	0,5 ml	100	C75VT1047

GLUCOSE

PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Potassium oxalate/Sodium fluoride - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, grey cap	0,5 ml	100	C75VT1048

COAGULATION

PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Sodium citrate 3,2%	13 x 75 mm, plastic, blue cap	1,8 ml	100	C75VT1013
	13 x 75 mm, DOUBLE WALL: inside tube PP, outside tube PET, blue cap	1,8 ml	100	C75VT1053



Sodium citrate 3,2%	13 x 75 mm, plastic, blue cap	2,7 ml	100	C75VT1014
	13 x 75 mm, plastic, blue cap	4,5 ml	100	C75VT1015



HEMATOLOGY



PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
EDTA K2	13 x 75 mm, plastic, purple cap	2 ml	100	C75VT1007
	13 x 75 mm, plastic, purple cap	3 ml	100	C75VT1008
	13 x 75 mm, plastic, purple cap	4 ml	100	C75VT1009
	13 x 100 mm, plastic, purple cap	6 ml	100	C75VT1051
EDTA K2 - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, purple cap	0,5 ml	100	C75VT1044

EDTA K3	13 x 75 mm, plastic, purple cap	2 ml	100	C75VT1010
	13 x 75 mm, plastic, purple cap	3 ml	100	C75VT1011
	13 x 75 mm, plastic, purple cap	4 ml	100	C75VT1012
	13 x 100 mm, plastic, purple cap	6 ml	100	C75VT1052
EDTA K3 - PEDIATRIC USE, micro tubes with integrated capillary	10 x 45 mm, plastic, purple cap	0,5 ml	100	C75VT1045

ESR


PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
ESR, reading on a rack with scale. FREE: at each 10 boxes x 100 tubes bought we grant for free a rack	8 x 120 mm, glass, black cap	1,6 ml	100	C75VT1016
	8 x 120 mm, glass, black cap	1,8 ml	100	C75VT1017
	8 x 120 mm, glass, black cap	2,0 ml	100	C75VT1018
ESR rack with scale, 10 positions	-	-	1	C66VT1019



OTHERS

PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Holders, for collection needles 18G, 20G, 21G, 22G and needles flashback or visio type 21G, 22G	-	-	1	EL310
	-	-	1	CZQ5






NEEDLES for VACUUM TUBES

PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Needles, sterile, non-pyrogenic, non-toxic, single use, safe and efficient, latex-free, allow multiple samples to be collected at a single puncture, sharp and smooth edges make painless penetration	1,20 mm x 38 mm, pink	18G	100	ELB18
	0,90 mm x 90 mm, yellow	20G	100	ELB20
	0,80 mm x 38 mm, green	21G	100	ELB21
	0,70 mm x 90 mm, black	22G	100	ELB22





NEEDLES for SYRINGES

PRODUCT NAME	DESCRIPTION	SIZE	SET	CODE
Needles, sterile, disposable, non-toxic, non-pyrogenic	1,27 mm x 38 mm, pink	18G	100	ELHN18
	0,90 mm x 38 mm, yellow	20G	100	ELHN20
	0,81 mm x 38 mm, green	21G	100	ELHN21
	0,71 mm x 32 mm, black	22G	100	ELHN22
	0,64 mm x 32 mm, blue	23G	100	ELHN23
	0,55 mm x 38 mm, purple	24G	100	ELHN24
	0,50 mm x 38 mm, orange	25G	100	ELHN25
	0,45 mm x 16 mm, brown	26G	100	ELHN26





🏢 **Tody Laboratories Int. SRL**
🏠 22 Vadul Moldovei, s1, OP 18
Bucharest, CP 014033, ROMANIA
📞 ISO 15189, 17043, 9001, 13485
🌐 www.todylaboratories.com

todylaboratories®
Responsible diagnostic. Since 1993.

www.todylaboratories.com

